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## Quellen

- Hoppe C, Poepel A, Elger CE. Epilepsy: Accuracy of patient seizure counts. Arch Neurol [Internet]. 2007;64(11):1595–9. Available from: <http://dx.doi.org/10.1001/archneur.64.11.1595>
- Beniczky S, Wiebe S, Jeppesen J, Tatum WO, Brazdil M, Wang Y, et al. Automated seizure detection using wearable devices: A clinical practice guideline of the International League Against Epilepsy and the International Federation of Clinical Neurophysiology. Clin Neurophysiol [Internet]. 2021;132(5):1173–84. Available from: <http://dx.doi.org/10.1016/j.clinph.2020.12.009>
- Fisher RS, Blum DE, DiVentura B, Vannest J, Hixson JD, Moss R, et al. Seizure diaries for clinical research and practice: limitations and future prospects. Epilepsy Behav [Internet]. 2012;24(3):304–10. Available from: <http://dx.doi.org/10.1016/j.yebeh.2012.04.128>
- Larsen PM, Beniczky S. Non-electroencephalogram-based seizure detection devices: State of the art and future perspectives. Epilepsy Behav [Internet]. 2023;148(109486):109486. Available from: <http://dx.doi.org/10.1016/j.yebeh.2023.109486>
- Arends J, Thijs RD, Gutter T, Ungureanu C, Cluitmans P, Van Dijk J, et al. Multimodal nocturnal seizure detection in a residential care setting: A long-term prospective trial. Neurology [Internet]. 2018;91(21):e2010–9. Available from: <http://dx.doi.org/10.1212/WNL.0000000000006545>
- van Westrhenen A, Lazeron RHC, van Dijk JP, Leijten FSS, Thijs RD, Dutch TeleEpilepsy Consortium. Multimodal nocturnal seizure detection in children with epilepsy: A prospective, multicenter, long-term, in-home trial. Epilepsia [Internet]. 2023;64(8):2137–52. Available from: <http://dx.doi.org/10.1111/epi.17654>
- Lazeron RHC, Thijs RD, Arends J, Gutter T, Cluitmans P, Van Dijk J, et al. Multimodal nocturnal seizure detection: Do we need to adapt algorithms for children? Epilepsia Open [Internet]. 2022;7(3):406–13. Available from: <http://dx.doi.org/10.1002/epi4.12618>
- Engelgeer A, van Westrhenen A, Thijs RD, Evers SMAA. An economic evaluation of the NightWatch for children with refractory epilepsy: Insight into the cost-effectiveness and cost-utility. Seizure [Internet]. 2022;101:156–61. Available from: <http://dx.doi.org/10.1016/j.seizure.2022.08.003>

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