
NightWatch +

- EN External connections for professional use
- NL Externe aansluitingen voor professioneel gebruik
- DE Externe Anschlüsse für den professionellen Einsatz
- FR Connexions externes pour une utilisation professionnelle
- IT Connessioni esterne per uso professionale
- ES Conexiones externas para el uso profesional

API v3

Legal Manufacturer

LivAssured BV

Schipholweg 103
2316 XC, Leiden
The Netherlands

M: info@nightwatchepilepsy.com

www.nightwatchepilepsy.com

© 2025 LivAssured, all rights reserved. No part of this manual may be reproduced in any form without the written permission of LivAssured.

LivAssured reserves the right to change or improve its products and any accompanying technical literature without specific notice of the changes or improvements.

External Connections Manual English Version 4.0

API version 3

Date of issue: 30th May 2025 (UTC/GMT +02:00 - Europe/Brussels)

Contents

1	General.....	2
1.1	About this manual.....	2
1.2	Explanation of symbols used in the manual.....	2
1.3	General warnings and cautions.....	2
1.4	Symbols on labels	3
2	External connections of NightWatch+	4
2.1	DIS relay connection.....	4
2.2	DAS API connection.....	4
3	DIS relay connection.....	6
4	DAS API connection.....	7
4.1	Configuration	7
4.2	Messages.....	8
4.3	Seizure alarm	9
4.4	Technical messages.....	9
4.5	Alarm station behavior	11
4.6	Testing the connection	11
4.7	Alarm delays	11
5	Alarm station parameters.....	11
6	Implementation	12
6.1	Client implementation	12
6.2	Server implementation.....	12
6.3	Infrastructure requirements.....	12
7	XML message examples.....	13
8	List of possible API message codes	16
9	Incident reporting.....	18
10	Contact information	18

1 General

1.1 About this manual

This external connections manual provides the information necessary to use NightWatch+ with the relay and API connections for professional use. It is complementary to the NightWatch+ User Manual. Please read the NightWatch+ User Manual before using NightWatch+. If any part of this manual is unclear, please contact support. The latest revision of the User Manual can be found at www.nightwatchchepilepsy.com

1.2 Explanation of symbols used in the manual

WARNING: Indicates a potentially hazardous situation which, if not avoided, could result in death or a serious injury.

CAUTION: Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the device.

PLEASE NOTE: A symbol used to emphasize information of which the user should be aware.

1.3 General warnings and cautions

PLEASE NOTE

- ▼ The Distributed Information System (DIS) relay connection and the Distributed Alarm System (DAS) API connection can only be used by professional users and are therefore disabled by default. The interfaces can only be enabled by the manufacturer. Contact LivAssured in this case.
- ▼ Connecting NightWatch+ to IT-networks including other equipment over the DIS or DAS connection could result in previously unidentified risks to patients, operators or third parties. The responsible organization must identify, analyze, evaluate and control their risks. Changes to the IT-network including network configuration, connecting additional items, disconnecting items, updating equipment and upgrading equipment could introduce new risks that require additional analysis.

CAUTION

- ▼ Only use the power adapters as supplied by LivAssured. Using a different charger or cable could damage the device and/or affect its performance.
- ▼ The parameters in the seizure alarm message are indicative and should not be used for medical treatment purposes.

WARNING

- ▼ The Distributed Information System (DIS) relay connection is a one-way forwarding of the alarm state of NightWatch+, intended to distribute information on the alarm state and should not be used to as a replacement for the visual and audible alarms. It cannot be relied upon as an alarming system as NightWatch+ cannot determine if information has been received by the DIS System or if the connection is lost.
- ▼ In case NightWatch+ is linked to a device that forwards alarms, be sure to regularly test if the link is working by following the instructions to manually triggering a NightWatch+ alarm and subsequently verify if the alarm is correctly transmitted and handled by the linked system. See "Section 4.6 Testing the connection".

1.4 Symbols on labels

The instruction manual must be read before use of the Device

The device contains electrical waste and must not be disposed of with normal household waste. The device must be disposed of in accordance with local law and the local code concerning electrical and electronic equipment.

Applied part, type BF

Manufactured in the Netherlands, manufacturing date

YYYY-MM

This symbol indicates a medical device

Serial number, production date followed by an identifier

Code to access the device monitoring data online

IP21

Classification of ingress protection by enclosure for alarm station, meaning: *Protected against access to hazardous parts with a finger and vertically falling drops of water or condensation*

IP22

Classification of ingress protection by enclosure for sensor, meaning: *Protected against access to hazardous parts with a finger and vertically falling drops of water when the enclosure is tilted at any angle up to 15° on either side of the vertical*

WARNING: Keep small parts away from children in the age range of 0 to 3 years. Spare and unassembled clips could pose a choking hazard when swallowed.

WARNING: Cords pose a strangulation hazard that may lead to death

CE Mark including notified body code

Temperature limits (minimum + maximum)

Pressure limits (minimum + maximum)

Humidity limits (minimum + maximum)

Keep dry

Keep away from heat

Direct current

PRO version marking

2 External connections of NightWatch+

NightWatch+ detects epileptic seizures and notifies caregivers using sound and light signals from the alarm station. The device also notifies caregivers in case of technical issues that make that the device cannot detect seizures. The device can also forward alarms or information to a connected (third party) information or alarming system.

The device offers two different ways to connect to third party information systems or alarm systems:

1. Using the RJ11 connector to connect to a Distributed Information System (DIS), using a dry-contact relay (**DIS relay connection**)
2. Using the RJ45 connector to connect to a Distributed Alarming System (DAS) using the API (**DAS API connection**).

These connections are enabled on PRO-versions of NightWatch+ devices. These devices can be identified by the PRO-version marking on the bottom of the alarm station:

2.1 DIS relay connection

A DIS connection is a very basic (digital) interface intended to distribute information to a network of connected devices which can forward this information to (remote) users (**DIS System**).

In NightWatch+, the DIS relay connection distributes information on the alarm state of NightWatch+ to the DIS system, which is either ON or OFF.

A DIS relay connection is a one-way information system. The RJ11 port of NightWatch+ is connected to an internal relay that is configured in a normally closed state. A device connected to the RJ11 port of the NightWatch+ can detect the opening of the relay, that represents an alarm state in the NightWatch+. There is no confirmation or guarantee that information distributed by NightWatch+ from the DIS relay connection to the network has been correctly received by the connected system.

The DIS relay can therefore only be used as a system that processes information on the alarm state of NightWatch+ in addition to the visual and audible alarms of the alarm station and cannot be used as a replacement for the visual and audible alarms of NightWatch+.

See "Section 3" for a technical description on how to connect NightWatch+ to a DIS system.

WARNING: The DIS relay connection is a one-way forwarding of the alarm state of NightWatch+, intended to distribute information on the alarm state and should not be used as a replacement for the visual and audible alarms. It cannot be relied upon as an alarming system as NightWatch+ cannot determine if information has been received by the DIS System or if the connection is lost.

2.2 DAS API connection

A DAS API connection is a digital network interface to distribute alarms or notifications to a device (**DAS system**). A DAS system can be programmed to forward messages received by the DAS API connection in different ways, or to different persons, depending on the urgency or severity of the event.

The DAS API connection serves as a two-way information system. A DAS system connected to the RJ45 port of NightWatch+ receives detailed information on the current state of the NightWatch+ and is notified on any events detected by NightWatch+. NightWatch+ receives a confirmation from the DAS System if a message has been received and will alarm if such a message could not be delivered.

The DAS System could serve as a replacement for the alarming function of the NightWatch+ alarm station, but only if the following conditions are met:

- The DAS is certified as a medical device class IIa or higher under (EU) 2017/745 Medical Device Regulation (MDR).
- The DAS should have a built-in alarm if no status updates have been received within a certain period of time.
- The DAS should be compliant with ISO 60601-1-8:2006 to forward alarming signals for high and medium priority alarms.
- The DAS should be compliant with ISO 80001-1:2021 risk management for IT-networks incorporating medical devices.

See "Section 4 DAS API connection" for a full description on how to connect the NightWatch+ to a DAS system.

PLEASE NOTE: The DIS relay connection and the DAS API connection can only be used by professional users and are therefore disabled by default. This functionality can only be enabled by the manufacturer. Contact LivAssured in this case.


WARNING: In case NightWatch+ is linked to a device that forwards alarms, be sure to regularly check if the link is working by following the instructions to manually trigger an NightWatch+ alarm and subsequently verify if the alarm is correctly transmitted and handled by the linked system. See Section 4 for a technical description on how to connect NightWatch+ to a DAS system.

WARNING: In case NightWatch+ is linked to a distributed messaging system that is not approved as a DAS according to these requirements, there is a risk that alarms do not reach the caregiver which may result in death.

3 DIS relay connection

The NightWatch+ alarm station has an RJ11 port, where the middle two of the pins serve as a dry-contact relay to forward information on the alarm state of the NightWatch+ alarm station. Below, the situation for a 4-pin connector is shown, with the middle two pins (2 & 3) are switched.

A 6-pin connector can also be used, the middle pins are in that case pins 3 & 4 that serve as the dry-contact relay.

The relay is closed when the NightWatch+ is not in an alarm state. This means that current can flow from pin 2 to pin 3. When NightWatch+ enters an alarm state, the NightWatch+ will open the relay for 2 seconds and then closes again.

The device will open the relay for 2 seconds when the following alarm states start:

Type	Description
Seizure Alarm	All seizure alarms will trigger the relay
Technical Alarm	The following technical alarms will trigger the relay: - Sensor connection lost - Sensor battery low - Sensor unable to track a heart rate for a while - Sensor switched off while being charged - Sensor accelerometer is defect - Alarm station lost power

The DIS relay connection is one-way information system, meaning that NightWatch+ cannot determine if information has been received by the DIS System. To ensure that the information is correctly received we advise to regularly provoke an alarm state on the NightWatch+ and check if the information is correctly received by the DIS system.

The DIS relay connection is only intended to be used by professional users that understand the limitations and risks of using a DIS system. It is not enabled by default. Contact LivAssured if you want to enable the DIS Relay Connection.

WARNING: The DIS relay connection is a one-way forwarding of the alarm state of NightWatch+, intended to distribute information on the alarm state and should not be used to as a replacement for the visual and audible alarms. It cannot be relied upon as an alarming system as NightWatch+ cannot determine if information has been received by the DIS System or if the connection is lost.

4 DAS API connection

The NightWatch+ alarm station can be connected to a network using the RJ45 ethernet connection.

Communication protocol

The communication setup uses a client-server model where the alarm station is the client initiating the connection. The remote server address needs to be configured in the alarm station. Contact LivAssured on how to do this.

The connection is set up by opening a TCP/IP Socket on the remote server. The client will normally only send relatively short messages and terminate the connection. However, it is possible that there are several clients trying to transmit simultaneously. The remote server must therefore listen on a specific port and support multiple parallel connections (multi-threading model).

4.1 Configuration

The Pro version of NightWatch+ includes a web-based interface. To access it, simply enter the alarm station's IP address into any browser connected to the same network.

The webpage will appear as shown below:

WARNING:

This web interface allows you to change the NightWatch settings which may affect the performance and function of the device.
The web interface should only be used by users that understand and accept the risks.
If in doubt, contact our helpdesk first.

ProVersion

Our version <currentversion> Remote Server has version <remoteversion>
Mac Address <##:##:##:##:##:##>
Sensor Version <moduleversion>

Administration Password

The text "ProVersion" is only visible on a pro-version of NightWatch+. If "ProVersion" is not visible, the DAS (and DIS relay) connection cannot be used. If "ProVersion" is not visible, but a pro-version of NightWatch+ was ordered, please contact LivAssured.

On the bottom of the alarm station on the label, there will be an 8-character access code.

Key : XXXXXXXX

On the web interface under "Administration Password", enter "DAS" followed by the 8-character code (DAS#####) to access the DAS configuration interface.

Network

The alarm station's network settings can be adjusted if required by the network infrastructure. In most cases, using DHCP is sufficient, and no changes are necessary.

Time

The alarm station needs access to a time server for event timestamps. If the alarm station has access to the internet, these settings do not require changing. If the network restrictions or policy require use of an internal time server, it can be configured here. You can check the system time at 'Board time in seconds'.

PRO Password

Use this option to change the device password for the web interface.

IpServer

This is the server to which API messages will be sent. To use the API, configure the server settings and ensure that 'Enable' is set to 'On'. The ClientAliveInterval determines how frequently this API message is sent, measured in seconds, and is set to 60 seconds by default." Finally, click "SET IP CONFIG" to save the settings.

4.2 Messages

The communication protocol defines the following steps:

- Client opens connection and sends announcement
- Server sends response
- Client sends one or more messages
- Server sends received message size
- Client confirms message size
- Client hangs up when done.

The format of the announcement:

```

struct announcement
{
    char    magic[16];        // Communication Identifier
    char    clientId[16];    // Client MAC address
    uint32_t version;        // API version
    uint32_t flags;          // Flags (supported items)
    uint64_t timestamp;      // Epoch Timestamp
}; // Total 48 bytes

```

The magic part of the announcement will be used to do a first rough check to see if the client is authorized to connect to the server.

Note: The Current Magic is "NightWatchAPI123"

- The alarm station MAC id will be used as the ClientID, in the format AA:BB:11:22:CC:DD.
- The current API version is 3
- There are no flags defined. Default is 0.
- Timestamp is in EPOCH, seconds since 1 jan 1970. (time(NULL))

The announcement must be followed by a message. The messages are defined as follows. Note that the message bodies will be sent in XML format.

```

struct message {
{
    uint32_t type;        // Message Type
    uint32_t size;        // Size of body, multiple of 8
    uint8_t body[];      // message body
}; // Total 8+size bytes

```

PLEASE NOTE: The Body of the message must be a multiple of 8 bytes.

The following different message body formats are defined:

Message	Frequency	Type	Remarks
Seizure Alarm	Directly	1	
Technical Message	Directly	2	
Client Alive	60 seconds	3	adjustable
(Reserved)	-	4	
Daily Status	When the sensor is placed on charger	5	

The message bodies contain the message in XML format. The details are given in the following subsections.

Message response

After receiving the announcement, message header and body, the Server must respond with 8 bytes:

- 4 bytes with "ACK\0", so string ACK terminated by NULL byte.
- 4 bytes with number of message bytes received. This should equal message size+8.

Client Alive

This message indicates that the client is still alive and up and running. Default settings are that it will be sent every 60 seconds. This setting can be changed in the PRO version or by LivAssured.

```
<clientalive>
<timestamp>20180810_223011</timestamp>
<epochtime>156234637</epochtime>
<macid>AA:BB:11:22:CC:DD</macid>
</clientalive>
```

Make sure to enable this function to ensure from a remote location that NightWatch+ is still online.

4.3 Seizure alarm

```
<seizurealarm>
<timestamp>20180810_223011</timestamp>
<epochtime>156234637</epochtime>
<macid> AA:BB:11:22:CC:DD</macid>
<code>8</code>
<desc>Motion3</desc>
<hr>69</hr>
</seizurealarm>
```

The Seizure Alarm message has the following format:

The macid provides the source MAC address of the alarm station generating the message.

The following seizure alarms are possible:

Description	Short	Code (hex)
Motion1	M1	2 (0x02)
Motion2	M2	4 (0x04)
Motion3	M3	8 (0x08)
PPRateIncrease	HI	16 (0x10)
HighPPRate	TC	32 (0x20)
LowPPRate	BC	128 (0x40)

These descriptions/codes refer to the different detection algorithms that are active in the sensor. PPM refers to the measured peaks-per-minute by the PPG sensor.

It is possible that multiple detection algorithms will trigger at the same time. In this case, the received code will be an OR operation of the individual detection algorithm. Example, if both Motion2 and PPM-Increase trigger the seizure detection, then the code $(4 | 16) = 20$ will be sent.

See section 8 for the full list of possible API codes.

CAUTION: The parameters in the seizure alarm messages are indicative and should not be used for medical treatment purposes.

4.4 Technical messages

```
<techalarm>
<timestamp>20180810_223011</timestamp>
<epochtime>156234637</epochtime>
<macid>AA:BB:11:22:CC:DD</macid>
<code>2048</code>
<desc>LowBattery</desc>
<hr>63</hr>
</techalarm>
```

The technical messages will be sent with the following XML body format:

The macid provides the source MAC address of the alarm station generating the message.

The following technical messages are possible.


Description	Short	Code	Explanation
Bad-Signal	BS	2	Sensor wearer is In-Rest and there is no heart rate signal, but there is motion, for >300 seconds or there is no Heart Rate signal and no motion for >120 seconds
All-Good	AG	65536	Sent as soon as the heart rate is detected confidentially (Standard detection mode)
Hw-Error	HW	8	Hardware Error, e.g. if the accelerometer is not functioning.
Start-Recording	SR	256	Sensor is removed from USB power
End-Recording	ER	512	Sensor is back on USB Power
Connection-Lost	CL	1024	Connection to the sensor is lost
Connection-Good	CG	2097152	The module connection has been restored, (after module lost).
Battery-Low	BL	2048	Battery is low, below 3.4 V
Module-Reboot	MR	4096	Module Reboot detected
Base-Reboot	BR	8192	Alarm Station was rebooted
Not-Charging	NC	16384	Sensor on charger, but sensor not powered on
(Reserved)		32768	(Reserved for future use)
In-Rest	IR	131072	Sent if patient is In-Rest state (Enhanced detection mode) meaning the sensor is positioned horizontally and no movement has been detected for >120 seconds
Power-Lost	PL	262144	Power of the alarm station is lost
Power-Good	PG	524288	Power of the alarm station is restored (after power lost)
Seizure-Ack	SA	1048576	Seizure alarm is acknowledged by pressing one of the buttons on the alarm station
IP server disconnected	IL	4194304	Connection with the third-party IP server lost. Typically, this would mean this message will not arrive, but it will try to send it regardless. This alarm state will also reset the alarm station audio to an audible level if the audio was turned down very low while connected to a DAS system.

See section 8 for the full list of possible API codes.

4.5 Alarm station behavior

Enabling DAS mode

When the alarm station has successfully connected with the DAS by sending an API message, it will switch to DAS mode. When DAS mode is enabled, the second blue LED will start to blink:

Volume setting

When DAS-mode is enabled, the alarm station volume can be lowered below the minimum volume set for regular operation. To do this, hold the volume down button () until the required volume is reached.

IP server lost alarm

When an API message cannot be delivered successfully to the DAS server, the alarm station will automatically exit DAS mode. When this happens, the volume will automatically be restored to the minimum volume set for regular operation. The alarm station will also give a technical alarm that the IP server connection has been lost, and a blinking technical alarm LED (orange) will be visible. When the DAS connection is restored, the alarm will stop and the volume will return to the level set in the original DAS setting, which could be lower than the minimum volume allowed in regular operation.

Disabling DAS mode

When the "IP server lost" alarm is active, press one of the volume buttons to disable DAS mode and stop the alarm.

DAS mode will automatically resume when a successful API message is received by the IP server.

4.6 Testing the connection

Technical alarm

When the API connection is implemented, a technical alarm message can be tested by triggering a Not-Charging (16384) message as follows:

- With the alarm station powered and sensor on the charger and switched ON, connect the API.
- With the sensor on the charger, switch the sensor off.
- Verify that the Not-Charging message with code 16384 is sent over the API

PLEASE NOTE: When the sensor is switched off, the alarm station will signal the "connection lost alarm". See the User Manual NightWatch+ for all the alarm station lights and sounds.

Epilepsy alarm

When the API connection is implemented, an epilepsy alarm can be tested by triggering as follows:

- With the alarm station powered and sensor on the charger and switched ON, connect the API.
- Take the sensor from the charger, the Start-Recording (256) message is sent over the API
- Shake the sensor rapidly back and forth until a seizure alarm is triggered. Verify that the alarm is Motion 3 (code 8).

4.7 Alarm delays

Alarms can have a delay of a maximum of 2 seconds from the onset of alarm condition to the point it leaves the output port.

5 Alarm station parameters

The following parameters are used on the alarm station.

- Remote Server IP Address
- Remote Server Port
- Encryption Enabled
- Type of messages included for seizure and technical messages

6 Implementation

Example code is available for both the client and the server.

6.1 Client implementation

The current client example implementation is done in pure C, using the standard sockets.

- Connection is opened to the server on given port
- The first announcement is sent with correct magic.
- Then XML message is sent encapsulated in the message struct (with size/type).

6.2 Server implementation

The current server example implementation is done in pure C++ with QT. It includes a simple graphical UI that displays the incoming messages. Main features:

- The server class listens for incoming connections (MwServer class using QTcpServer).
- If connection is coming in, a new thread is started that will handle the connection (NwSrvThread class using QThread).
- Each server thread will extract the incoming XML message and forward them to top level application (using QT's Signal/Slot mechanism). NwSrvThread will hang-up (and delete itself) if no new data is received withing 10 seconds.

Note that third parties that wish to implement their own server can use the QT NwSrv as an example, and the client implementation as a test program.

6.3 Infrastructure requirements

The server can be configured with any OS and approved alarming software if it can communicate with the xml message structure over the TCP/IP protocol. Multiple NightWatch+ systems can be connected to the same network. The maximum message size is 242 bytes, so the network bandwidth reserved for the connected NightWatch+ systems will have to be at least 242 bytes times the number of systems.

7 XML message examples

The following is a list of XML Message body examples:

```
<clientalive>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
</clientalive>
```

```
<seizurealarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>2</code>
  <desc>Motion1</desc>
</seizurealarm>
```

```
<seizurealarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>4</code>
  <desc>Motion2</desc>
</seizurealarm>
```

```
<seizurealarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>8</code>
  <desc>Motion3</desc>
</seizurealarm>
```

```
<seizurealarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>16</code>
  <desc>PPRateIncrease</desc>
</seizurealarm>
```

```
<seizurealarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>32</code>
  <desc>HighPPRate</desc>
</seizurealarm>
```

```
<seizurealarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>128</code>
  <desc>LowPPRate</desc>
</seizurealarm>
```

```
<techalarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>2</code>
  <desc>Bad-Signal</desc>
</techalarm>
```

```
<techalarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>65536</code>
  <desc>All-Good</desc>
</techalarm>
```

```
<techalarm>
  <timestamp>20201109_220911</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>8</code>
  <desc>Hw-Error</desc>
</techalarm>
```

```
<techalarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>256</code>
  <desc>Start-Recording</desc>
</techalarm>
```

```
<techalarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>512</code>
  <desc>End-Recording</desc>
</techalarm>
```

```
<techalarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>1024</code>
  <desc>Connection-Lost</desc>
</techalarm>
```

```
<techalarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>2097152</code>
  <desc>Connection-Good</desc>
</techalarm>
```

```
<techalarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>2048</code>
  <desc>Battery-Low</desc>
</techalarm>
```

```
<techalarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>4096</code>
  <desc>Module-Reboot</desc>
</techalarm>
```

```
<techalarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>8192</code>
  <desc>Base-Reboot</desc>
</techalarm>
```

```
<techalarm>
  <timestamp>20180810_223011</timestamp>
  <epochtime>156234637</epochtime>
  <macid>AA:BB:11:22:CC:DD</macid>
  <code>16384</code>
  <desc>Not-Charging</desc>
</techalarm>
```

```
<techalarm>  
  <timestamp>20180810_223011</timestamp>  
  <epochtime>156234637</epochtime>  
  <macid>AA:BB:11:22:CC:DD</macid>  
  <code>131072</code>  
  <desc>In-Rest</desc>  
</techalarm>
```

```
<techalarm>  
  <timestamp>20180810_223011</timestamp>  
  <epochtime>156234637</epochtime>  
  <macid>AA:BB:11:22:CC:DD</macid>  
  <code>262144</code>  
  <desc>Power-Lost</desc>  
</techalarm>
```

```
<techalarm>  
  <timestamp>20180810_223011</timestamp>  
  <epochtime>156234637</epochtime>  
  <macid>AA:BB:11:22:CC:DD</macid>  
  <code>524288</code>  
  <desc>Power-Good</desc>  
</techalarm>
```

```
<techalarm>  
  <timestamp>20180810_223011</timestamp>  
  <epochtime>156234637</epochtime>  
  <macid>AA:BB:11:22:CC:DD</macid>  
  <code>1048576</code>  
  <desc>Seizure-Ack</desc>  
</techalarm>
```

8 List of possible API message codes

Type	Code	Suggested description	Breakdown of the code
Seizure	2	Seizure alarm: Motion	Motion1
Seizure	3	Seizure alarm: Motion	Motion1 + UprightPosition
Seizure	4	Seizure alarm: Motion	Motion2
Seizure	5	Seizure alarm: Motion	Motion2 + UprightPosition
Seizure	6	Seizure alarm: Motion	Motion1 + Motion2
Seizure	7	Seizure alarm: Motion	Motion1 + Motion2 + UprightPosition
Seizure	8	Seizure alarm: Motion	Motion3
Seizure	10	Seizure alarm: Motion	Motion1 + Motion3
Seizure	12	Seizure alarm: Motion	Motion2 + Motion3
Seizure	14	Seizure alarm: Motion	Motion1 + Motion2 + Motion3
Seizure	16	Seizure alarm: High PPM	PPMIncrease
Seizure	18	Seizure alarm: High PPM	PPMIncrease + Motion1
Seizure	20	Seizure alarm: High PPM	PPMIncrease + Motion2
Seizure	22	Seizure alarm: High PPM	PPMIncrease + Motion1 + Motion2
Seizure	24	Seizure alarm: High PPM	PPMIncrease + Motion3
Seizure	26	Seizure alarm: High PPM	PPMIncrease + Motion1 + Motion3
Seizure	28	Seizure alarm: High PPM	PPMIncrease + Motion2 + Motion3

Type	Code	Suggested description	Breakdown of the code
Seizure	30	Seizure alarm: High PPM	PPMIncrease + Motion1 + Motion2 + Motion3
Seizure	32	Seizure alarm: High PPM	HighPPM
Seizure	34	Seizure alarm: High PPM	HighPPM + Motion1
Seizure	36	Seizure alarm: High PPM	HighPPM + Motion2
Seizure	38	Seizure alarm: High PPM	HighPPM + Motion1 + Motion2
Seizure	40	Seizure alarm: High PPM	HighPPM + Motion3
Seizure	42	Seizure alarm: High PPM	HighPPM + Motion1 + Motion3
Seizure	44	Seizure alarm: High PPM	HighPPM + Motion2 + Motion3
Seizure	46	Seizure alarm: High PPM	HighPPM + Motion1 + Motion2 + Motion3
Seizure	48	Seizure alarm: High PPM	PPMIncrease + HighPPM
Seizure	50	Seizure alarm: High PPM	PPMIncrease + HighPPM + Motion1
Seizure	52	Seizure alarm: High PPM	PPMIncrease + HighPPM + Motion2
Seizure	54	Seizure alarm: High PPM	PPMIncrease + HighPPM + Motion1 + Motion2
Seizure	56	Seizure alarm: High PPM	PPMIncrease + HighPPM + Motion3
Seizure	58	Seizure alarm: High PPM	PPMIncrease + HighPPM + Motion1 + Motion3
Seizure	60	Seizure alarm: High PPM	PPMIncrease + HighPPM + Motion2 + Motion3
Seizure	62	Seizure alarm: High PPM	PPMIncrease + HighPPM + Motion1 + Motion2 + Motion3

Type	Code	Suggested description	Breakdown of the code
Seizure	128	Seizure alarm: Low PPM	LowPPM
Seizure	130	Seizure alarm: Low PPM	LowPPM + Motion1
Seizure	132	Seizure alarm: Low PPM	LowPPM + Motion2
Seizure	134	Seizure alarm: Low PPM	LowPPM + Motion1 + Motion2
Seizure	136	Seizure alarm: Low PPM	LowPPM + Motion3
Seizure	138	Seizure alarm: Low PPM	LowPPM + Motion1 + Motion3
Seizure	140	Seizure alarm: Low PPM	LowPPM + Motion2 + Motion3
Seizure	142	Seizure alarm: Low PPM	LowPPM + Motion1 + Motion2 + Motion3
Tech	2	Technical issue: Unable to detect heart rate	Bad-Signal
Tech	8	Technical issue: Hardware error	Hw-Error
Tech	256	NightWatch recording started	Start-Recording
Tech	512	NightWatch recording stopped	End-Recording
Tech	1024	Technical issue: Sensor connection lost	Connection-Lost
Tech	2048	Technical issue: Sensor battery low	Battery-Low

Type	Code	Suggested description	Breakdown of the code
Tech	4096	Sensor rebooted	Module-Reboot
Tech	8192	Base rebooted	Base-Reboot
Tech	16384	Technical issue: Sensor switched off	Not-Charging
Tech	32768	-	(Reserved)
Tech	65536	Wearer not in rest – low sensitivity	All-Good
Tech	131072	Wearer in rest – high sensitivity	In-Rest
Tech	262144	Power Lost	Power-Lost
Tech	524288	Power Good	Power-Good
Tech	1048576	Seizure Acknowledged	Seizure-Ack
Tech	2097152	Module connection restored	Connection-Good
Tech	4194304	Connection with the third party IP server lost	IP server disconnected

9 Incident reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the country in which the user is established.

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Support: +3185 0601252

Email: info@nightwatchepilepsy.com

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ProVersion

Our version <currentversion> Remote Server has version <remoteversion>
Mac Address <###:###:###:###:###:###>
Sensor Version <moduleversion>

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The screenshot shows the NightWatch+ web interface with the following sections:

- Sensors:** HeartRate Sensor (Port 7) ONLINE Battery: 4.1 V Sensor Version: (Software version)
- Control:**
 - Mac Address: (mac address)
 - Board time in seconds: (System time) 0:00
 - Network:**
 - DHCP: ☒ ON
 - IP Address:
 - Netmask:
 - Gateway:
 - Dns Server:
 - Hostname:
 - [SET NETWORK](#)
 - Time:**
 - Server: ipennightwatch.nl
 - Pool Interval: 600
 - [SET TIME SERVER](#)
 - Pro:**
 - Password:
 - [SET PRO PASSWORD](#)
 - IpServer:**
 - Enable: ☐ OFF
 - IP Address:
 - Port:
 - ClientAliveInterval: 60
 - [SET IP CONFIG](#)
- Set Config Defaults (will reboot base):** [>](#)
- Reboot Alarm Station:** [>](#)
- Log Out:** [>](#)

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The screenshot shows the 'Sensors' configuration window for a 'HeartRate Sensor (Port 7)'. The status is 'ONLINE' with a battery level of 4.1 V. The interface is divided into several sections: 'Control' with a 'Board time in seconds' display; 'Network' with fields for IP Address, Netmask, Gateway, Dns Server, and Hostname, and a 'DHCP' toggle set to 'ON'; 'Time' with fields for Server, Poll Interval, and a 'Set Time Server' button; 'Pro' with a 'Password' field and a 'Set Pro Password' button; and 'IpServer' with an 'Enable' toggle set to 'OFF' and fields for IP Address, Port, and ClientAliveInterval. At the bottom, there are links for 'Set Config Defaults (will reboot base)', 'Reboot Alarm Station', and 'Log Out'.

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