NightWatch +

EN User manual

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M: info@nightwatchepilepsy.com

www.nightwatchepilepsy.com

© 2025 LivAssured, all rights reserved. No part of this manual may be reproduced in any form without the written permission of LivAssured.

LivAssured reserves the right to change or improve its products and any accompanying technical literature without specific notice of the changes or improvements.

User Manual English Version 7.0

Date of issue: 2nd April 2025 (UTC/GMT +02:00 - Europe/Brussels)

User Manual

EN

Thank you for choosing the nocturnal epileptic seizure detection and monitoring device,

NightWatch +

We at LivAssured | NightWatch understand that caring for persons suffering from epileptic seizures is highly demanding and stressful.

Seizures can be scary, as some of them can result in injuries or even sudden unexpected death in epilepsy (SUDEP) especially when a person suffering from seizures is unattended at night. NightWatch+ will warn you for the most dangerous seizures during sleep and it has been proven to reduce stress for the caregiver.

NightWatch was invented because multiple neurologists from the Dutch Academic Centre for Epileptology at Kempenhaeghe and Epilepsy Centre SEIN saw that there was a need for a reliable epileptic seizure detection device. These neurologists started a cooperation, a so-called consortium, of Dutch neurologists from

Kempenhaeghe and SEIN, multiple universities, and patient organizations. This consortium invented, developed and clinically validated the first version of the NightWatch. LivAssured | NightWatch was founded and further developed the NightWatch in cooperation with the consortium. This resulted in the high-quality and reliable device in front of you today.

LivAssured | NightWatch has the mission to improve the lives of people living with epilepsy. We do this by developing products like NightWatch+ that support the daily care for people with epilepsy and which enable research into new and better treatments.

In this manual, we explain how to get started and use your NightWatch+.
We stand beside you with our product, as well as with help, advice and tips.

Contents

1	General.....	4	6.2	Information signals	19
1.1	About this Manual.....	4	6.3	Responding to alarms.....	19
1.2	Explanation of symbols used in the manual.....	4	7	Seizure monitoring.....	22
1.3	General warnings and cautions.....	4	7.1	The NightWatch Portal.....	22
1.4	Symbols on labels	5	7.2	How to connect to the NightWatch Portal.....	22
2	NightWatch+	6	7.3	Does the NightWatch Portal respect my privacy?.....	22
2.1	Intended use.....	6	7.4	Connecting NightWatch+ to alarm systems.....	22
2.2	Medical classification	6	8	Specifications.....	23
2.3	Type of seizures NightWatch+ detects.....	6	8.1	Technical Specifications	23
2.4	Operation mode.....	6	8.2	Environmental conditions.....	23
2.5	User profile	7	8.3	Electromagnetic environment conditions.....	24
2.6	Possible side effects.....	7	8.4	Regulatory and Compliance.....	27
2.7	Clinical performance & benefits	7	8.5	Cybersecurity.....	27
3	NightWatch+ package contents	9	9	Incident reporting	28
4	The different parts and how they work	10	10	Service life and guarantee	28
4.1	Sensor with elastic strap.....	10	11	Disposal.....	28
4.2	Comfort Patch.....	10	12	Contact information	28
4.3	Alarm station.....	11			
5	Using your NightWatch+	12			
5.1	Preparing for first use	12			
5.2	Daily use	15			
5.3	Managing alarm sounds.....	16			
5.4	Charging the sensor.....	17			
5.5	Maintenance & cleaning	17			
5.6	Reuse.....	17			
5.7	Transport or storage.....	17			
6	Signals from the sensor and the alarm station	18			
6.1	Alarm signals	18			

1 General

1.1 About this Manual

This manual provides the information necessary to use NightWatch+ safely and effectively and is available in multiple languages. Please read the manual before using NightWatch+. If any part of this manual is unclear, please contact support. See Section 12 for contact details. The latest version of the User Manual is available at www.nightwatchepilepsy.com

1.2 Explanation of symbols used in the manual

WARNING: Indicates a potentially hazardous situation which, if not avoided, could result in death or a serious injury.

CAUTION: Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the device.

PLEASE NOTE: A symbol used to emphasize information of which the user should be aware.

1.3 General warnings and cautions

CAUTION

- NEVER diagnose or treat yourself based on the readings of NightWatch+. ALWAYS consult with your physician.
- Do not accept and use the device if there are signs of piercing, manipulation, water damage or any other damage to the device, its packaging or label. Contact your supplier for help.
- Only use the power adapters as supplied by LivAssured. Using a different charger or cable could damage the device and/or affect its performance.

WARNING

- This product does not guarantee that 100% of the epileptic seizures will be detected, therefore there is a possibility that the caregiver is not warned of an epileptic seizure while using NightWatch+
- No modification of this equipment is allowed. Modifications to the device could lead to the hardware, algorithm, connectivity or communication to fail.
- Use of transducers and cables other than those specified or provided by LivAssured could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of NightWatch+, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Keep small parts away from children in the age range of 0 to 3 years. Spare and unassembled clips could pose a choking hazard when swallowed.
- Cords pose a strangulation hazard that may lead to death.

1.4 Symbols on labels

The instruction manual must be read before use of the Device

The device contains electrical waste and must not be disposed of with normal household waste. The device must be disposed of in accordance with local law and the local code concerning electrical and electronic equipment.

Applied part, type BF

Manufactured in the Netherlands, manufacturing date

YYYY-MM

This symbol indicates a medical device

Serial number, production date followed by an identifier

Code to access the device monitoring data online

IP21

Classification of ingress protection by enclosure for alarm station, meaning: *Protected against access to hazardous parts with a finger and vertically falling drops of water or condensation*

IP22

Classification of ingress protection by enclosure for sensor, meaning: *Protected against access to hazardous parts with a finger and vertically falling drops of water when the enclosure is tilted at any angle up to 15° on either side of the vertical*

WARNING: Keep small parts away from children in the age range of 0 to 3 years. Spare and unassembled clips could pose a choking hazard when swallowed.

WARNING: Cords pose a strangulation hazard that may lead to death

CE Mark including notified body code

Temperature limits (minimum + maximum)

Pressure limits (minimum + maximum)

Humidity limits (minimum + maximum)

Keep dry

Keep away from heat

Direct current

2 NightWatch+

2.1 Intended use

NightWatch+ is intended to notify a caregiver of the occurrence of a patient's Nocturnal Epileptic Motor Seizures(*) and monitor seizure frequency over time.

(*)Nocturnal Epileptic Motor Seizures, being the following seizure types¹:

- ▼ Tonic-clonic
- ▼ Tonic (if clustered or prolonged)
- ▼ Myoclonic (if clustered)
- ▼ Hyperkinetic

2.2 Medical classification

NightWatch+ is a class IIa device in accordance with the rules of the EU regulation concerning medical devices, (EU)2017/745, dated 5 April 2017 (the "Medical Device Regulation").

2.3 Type of seizures NightWatch+ detects

Seizures taking place during the night/while sleeping are termed nocturnal seizures. Motor seizures are any type of seizures involving muscle activity, such as sudden stiffness or tension in the muscles of the arms, legs, or trunk during a tonic seizure, or leg pedaling movements during a hyperkinetic seizure.

The heart rate is controlled by the autonomic nervous system. Epileptic seizures affect this system in complex ways. Increases in heart rate are common during motor seizures but sudden decreases in heart rate are also observed. NightWatch uses these changes in heart rate to detect seizures.

NightWatch+ is intended to notify for the most dangerous nocturnal motor seizures associated with a risk of Sudden Unexpected Death in Epilepsy (SUDEP) or injury², which are the seizure types described in the intended use.

These seizures are detected by combining and analyzing heart rate data, measured with a photoplethysmography (PPG) sensor, and movement data, measured with an accelerometer (ACC).

2.4 Operation mode

NightWatch+ is a wearable device consisting of a wireless sensor and an alarm station. The sensor is worn during sleep on the biceps of the upper arm. It includes a PPG (photoplethysmography) sensor to track the heart rate, an ACC (accelerometry) movement sensor, a microprocessor that processes the data from the sensors using a detection algorithm and a battery.

NightWatch+ does not provide direct monitoring of the tracked heart rate or movement data. NightWatch+ is not a heart rate monitor.

The detection algorithm detects if the sensor readings match preprogrammed parameters that are associated with Nocturnal Epileptic Motor Seizures. If a match is detected, the epilepsy alarm is triggered and transmitted to the alarm station. It is neither possible nor necessary to modify the algorithms or adjust any thresholds to influence the device's performance.

The sensor and alarm station communicate using a wireless DECT protocol. When an epilepsy alarm is transmitted from the sensor to the alarm station, the alarm station notifies caregivers with an audible alarm and a blinking red LED light. A caregiver can then go to the person with epilepsy and, if necessary, provide assistance according to their physician's instructions.

The alarm station also warns caregivers with alarm sounds and orange blinking LED lights in case the system is unable to detect seizures for technical reasons. Possible reasons include a depleted battery, a lost connection between the sensor and alarm station (out of range), or the sensor being unable to track a PPG signal or movement data needed for seizure detection.

CAUTION: NEVER diagnose or treat yourself based on the readings of NightWatch+. ALWAYS consult with your physician.

PLEASE NOTE: Do not use this product in an environment where DECT signals may be blocked or interfered with by environmental properties or other equipment.

¹ Nomenclature is based on the classification of epileptic seizures by the International League Against Epilepsy (ILAE): S. Fisher, J. H. Cross, C. D'Souza, J. A. French, S. R. Haut, N. Higurashi, et al. "Instruction manual for the ILAE 2017 operational classification of seizure types," *Epilepsia*, vol. 58, no. 4, pp. 531-542, 2017

² J. Arends, R. D. Thijs, T. Gutter, C. Ungureanu, P. Cluitmans, J. Van Dijk, J. van Andel, F. Tan, A. de Weerd, B. Vledder, W. Hofstra, R. Lazeron, G. van Thiel, K. C. B. Roes and F. Leijten, "Multimodal nocturnal seizure detection in a residential care setting: A long-term prospective trial," *Neurology*, vol. 91, no. 21, pp. e2010-e2019, 2018

2.5 User profile

Intended users

Users of NightWatch+ are people diagnosed with epilepsy aged 4 years and older, having Nocturnal Epileptic Motor Seizures and caregivers thereof.

PLEASE NOTE: Patients using NightWatch+ can act as operator of the device during setup and usage, but someone has to be able to respond to alarms when a seizure is detected.

Contra-indications

NightWatch+ has no absolute contraindications that make its use completely inadvisable. However, certain risk factors do not necessarily exclude a user but require increased attention and may make the use of NightWatch+ inadvisable.

- ▼ User with damaged skin on the upper arm(s)

CAUTION: NightWatch+ is found to be biological safe to be worn on intact skin. Wearing the device on damaged skin could cause (further) irritation or damage of the contact area.

- ▼ User with cardiac arrhythmia

CAUTION: NightWatch+ detects seizures by measuring, among other parameters, sudden changes in heart rate. Cardiac arrhythmia could influence the performance of NightWatch+.

- ▼ User aged younger than 4 years

CAUTION: NightWatch+ is not intended for children younger than 4 years old. For children younger than 4 years old, clinical data evaluating the possible risks and benefits of using the device is lacking.

Intended use environment

NightWatch+ is intended to be used at home or at residential care facilities. NightWatch+ is not intended to be used in intensive care environments.

2.6 Possible side effects

For the previous generation of NightWatch+, the NightWatch Original (which consists of the same materials), it was observed that less than 3.5% of the users developed skin irritation when wearing the sensor without the comfort patch. See Chapter '4.2 Comfort Patch'.

2.7 Clinical performance & benefits

Epileptic seizures can result in injuries, status epilepticus, and even Sudden Unexpected Death in Epilepsy (SUDEP).

Scientific studies have shown that nocturnal tonic-clonic seizures are the most dangerous type that can lead to injuries or hospitalization. NightWatch+ uses advanced technology and algorithms, previously validated in multiple scientific studies, to detect these seizures which enables that assistance can be provided^{2,3,4}.

Seizure detection sensitivity

NightWatch+ has a seizure detection sensitivity for Nocturnal Epileptic Motor Seizures and for tonic-clonic seizures in adults and children (ages 4-18) which are respectively:

- ▼ Adults: 86% median (95% CI: 71%-93%) for Nocturnal Epileptic Motor seizures and 96% median (95% CI: 80%-100%) for tonic-clonic seizures
- ▼ Children: 100% median (95% CI: 87%-100%) for Nocturnal Epileptic Motor seizures and 100% median (95% CI*: 100%-100%) for tonic-clonic seizures

Using NightWatch+ leads to reduced stress for caregivers of people with epilepsy.

WARNING: This product does not guarantee that 100% of the epileptic seizures will be detected, therefore there is a possibility that the caregiver is not warned of an epileptic seizure while using NightWatch+.

³ R. Lazeron, R. Thijs, J. Arends, T. Gutter, P. Cluitmans, J. Van Dijk, F. Tan, W. Hofstra, C. Donjacour and F. Leijten, "Multimodal nocturnal seizure detection: Do we need to adapt algorithms for children?," Epilepsia, vol. 7, no. 3, pp. 406-413, 2022

⁴ A. van Westrhenen, R. H. C. Lazeron, J. P. van Dijk, F. S. S. Leijten and R. D. Thijs, "Multimodal nocturnal seizure detection in children with epilepsy: A prospective, multicenter, long-term, in-home trial," Epilepsia, vol. 64, pp. 2137-2152, 2023.

False alarm rate

NightWatch+ may sometimes trigger an epilepsy alarm when no epileptic seizure is occurring or it triggers an alarm for a minor seizure, which is a seizure of very short duration. When an epilepsy alarm was triggered by a minor seizure, it may not be visibly apparent that a seizure occurred. These occurrences are called false alarms.

NightWatch+ has a false alarm rate for Nocturnal Epileptic Motor Seizure detection in adults and children (ages 4-18), which are respectively:

- ▼ Adults: 0.03 per hour (median)
- ▼ Children: 0.04 per hour (median)

This means that most users will experience fewer than 1 false alarm every 3 to 4 nights.

CAUTION: This device may also give seizure alarms if no seizure is taking place which could cause stress for the caregiver.

PLEASE NOTE: Please contact LivAssured's service department when you have many false alarms to try to help you resolve this. See section 12 for contact details.

Clinical Testing details

NightWatch was invented, developed, and clinically validated by a consortium of Dutch neurologists from the Academic Center for Epileptology at Kempenhaeghe and epilepsy center SEIN, multiple universities, and patient organizations. The device was tested in scientific studies both in institutions and in home environments.

	Neurology ^a	Epilepsia ^a
	Arends et al. 2018 ²	Westrhenen et al 2023 ⁴
Patients	28	53
Age	15-67	4-16
Location	Institution	Home
Nights	1826	2310
Seizures	809	552
Sensitivity tonic-clonic seizures (median)	96% (95% CI*: 80%-100%)	100% (95% CI*: 100%-100%)
Sensitivity All seizure types (median)	86% (95% CI*: 77%-93%)	100% (95% CI*: 87%-100%)
False alarm rate/ hour	0.03	0.04

*: 95% Confidence Interval means that if the same population were to be sampled on multiple occasions, for 95 percent of the cases the median sensitivity will fall within the range.

3 NightWatch+ package contents

- 1 Alarm station
 - ▼ Includes built-in backup battery: Lithium-ion battery 3.7V, 450mAh, not replaceable
- 2 The FRIWO NEO006.0-I-X-05, 5VDC/1.4A power adapter with barrel jack for the alarm station (black)
- 3 Sensor clips (1 already in the sensor, 3 additional including 2 spare clips)
- 4 Comfort patches (3x)
- 5 Travel case
- 6 Sensor
 - ▼ Operates on built-in battery: Lithium-ion battery 3.7V, 450mAh, not replaceable
- 7 Sensor elastic strap (1 meter)
- 8 The FRIWO NEO006.0-I-X-05, 5VDC/1.4A power adapter with USB-C connector for the sensor (black)
- 9 Manual
- 10 Quick guides



4 The different parts and how they work

4.1 Sensor with elastic strap

The sensor is worn on the upper arm and secured with an elastic strap. The optimal position for the sensor is around the upper arm on the front of the biceps, not on the side. This positioning helps prevent the wearer from lying on the sensor when turning onto their side, as doing so could interfere with heart rate tracking. The orientation of the indicator light, whether facing up or down, does not affect functionality.

- A = ON/OFF switch
- B = PPG sensor & LEDs
- C = Charging port
- D = Indicator light

ON/OFF switch

The dark grey underside of the sensor contains an ON/OFF switch (A). The position can be switched by sliding it to the side using a pointy object (pen).

The ON/OFF switch is usually only used for the first activation.

You can use the ON/OFF switch to turn OFF the sensor completely when the sensor cannot be charged and is not being used. This is to prevent the battery from depleting.

Heart rate and movement sensor

The sensor continuously tracks the wearer's heart rate and movement to detect epileptic seizures. It uses two green LEDs (B) on the dark grey bottom of the sensor to track the heart rate through Photoplethysmography (PPG). One green heart rate LED can shut itself off when there is too much environmental light. Both LEDs should turn on in the first seconds after you remove the sensor from the charger.

Charging port

The charging port (C) is used to charge the sensor with the supplied USB-C power adapter.


Indicator light

The notch on the white side of the sensor contains an indicator light (D). The brightness of this indicator is set for use in a darkened room and may consequently be more difficult to see in daylight / a brightly lit area. The indicator light is used for several signals. Read about these signals in Chapter 6.

4.2 Comfort Patch

The highest-grade materials were selected when NightWatch+ was designed. The sensor and strap have been subjected to biocompatibility analysis and have been found biologically safe for its use according to the ISO 10993 standard. It is therefore unlikely that you'll develop an allergic reaction from NightWatch+.

However, it was observed for the previous generation of NightWatch+, the NightWatch Original (which consists of the same materials) that less than 3.5% of the users developed skin irritation when wearing the sensor without the comfort patch. We therefore recommend applying the comfort patch to offer the highest comfort while wearing the sensor.

Please make sure to replace the comfort patch regularly when it becomes loose, damaged or dirty. We recommend replacing it after one month. More comfort patches can be ordered via the website or by contacting LivAssured. See section 12 for contact details.

CAUTION: Please make sure to replace the comfort patch regularly when it becomes loose, damaged or dirty. Due to the buildup of sweat and bacteria on the comfort patch, skin irritation is more likely to occur.

CAUTION: The sensor of NightWatch+ should only be worn on intact skin. Do not continue using the device on the same skin location in case the skin becomes red, itchy or if any pain is felt and place a comfort patch on the NightWatch+ sensor. Please contact LivAssured in this case. See section 12 for contact details.

Applying the Comfort Patch


Remove the backing paper from the bottom of the patch and adhere it to the dark side of the NightWatch+ sensor. Use the middle hole and the opening for the switch as position guides, as shown in the figure below.

- ▼ Make sure that the charging port remains free for charging.
- ▼ Make sure that the ON/OFF switch of the sensor remains free.
- ▼ Make sure the patch does not cover the black sensor area in the middle.

4.3 Alarm station

The square white box is the alarm station. It emits both light and sound signals. The alarm station should be placed indoors near a power outlet.


A = Loudspeaker

B = LED indicator lights

C = Volume controls

🔊 Decreases the volume

🔊 Increases the volume

D = Network connection (RJ45 port)

E = RJ11 port (disabled, professional use only)

F = Power supply

Alarm station icons

- 🔌 LED 1 (blue) – Alarm station power status
- 📶 LED 2 (blue) – Network status
- 🔌 LED 3 (blue) – Sensor charging status
- 🔊 ● 🛑 LED 4 (green) – Audio paused state / Operation mode
- 🔔 LED 5 (orange) – Technical alarm
- 🚨 LED 6 (red) – Epilepsy alarm

5 Using your NightWatch+

Scan this QR code to watch a video online explaining the setup of the NightWatch and how to start using it.

5.1 Preparing for first use

Step 1: Adjust the sensor strap

You will need:

- ▼ Elastic strap, with clip inserted in sensor
- ▼ Additional clip
- ▼ Sensor
- ▼ Ballpoint pen/pencil
- ▼ Sharp (fabric) scissors

1

The elastic strap should be adjusted to the size of the upper arm of the wearer. One side of the elastic strap has already been fitted into the sensor.

2


Use the remaining end of the elastic strap to measure the circumference of the wearer's arm above the biceps in a bent position. Do not stretch the strap.

3

Mark the elastic strap where it overlaps the second sensor clip entrance and cut the elastic strap at the marked point.

4


The sensor must not fit too tight but should fit snugly against the skin.

5


Cut the strap on the marked location.

6

Attach a clip to the remaining end of the strap by pushing one corner of the elastic strap through the flat side of the clip. Pull the strap through the clip until it sticks out at a length of around one centimeter.


7

Push the clip into the remaining slot of the sensor and make sure that the end of the strap is protruding from the bottom of the sensor.

8

Place the sensor around the wearer's upper arm on the front of the biceps, not to the side. This ensures that the wearer is unlikely to lie on the sensor when turning onto his/her side. Lying on the sensor could disrupt tracking the heart rate. It does not matter if the indicator light is facing up or down.

PLEASE NOTE: In order to achieve as accurate detection as possible of any epileptic seizures, the sensor should be worn on the upper arm (biceps) and directly on the skin. If the wearer is wearing a pajama top or other long-sleeved garment, please ensure the sleeve is wide enough to wear the sensor underneath or else seizure detection is not possible.

Check the tightness of the strap. The correct fit allows for just about one finger underneath the strap.

9

If the elastic strap feels too loose:

Pull the elastic strap further through the clips. Try again.

If the elastic strap feels too tight:

- Using a fingernail, press the clips from below to remove them from the sensor.
- Adjust the clips to make the space between the clips on the elastic strap longer.
- Place the clips back into the slots and try again.

Adjust the strap until it is both comfortable for the wearer and fits snugly around the arm. Finally, cut off the excess ends of the strap with scissors to ensure it does not cover the green PPG LEDs and sensor.

There is no cause for concern if the sensor leaves an imprint on the arm after a sleep cycle, as long as this imprint fades by itself within hours.

PLEASE NOTE: Please ensure the sensor has been adjusted to the correct size for the patient before use. If it fits too tight or too loose, seizure detection may be less accurate.


Please ensure that the protruding ends of the elastic strap do not cover the green LEDs of the PPG sensor or seizure detection may be less accurate.

Step 2: First activation of the sensor

- 1 Take a pointed object or a pen/pencil
- 2 The dark underside of the sensor houses an ON/OFF switch. Use a pencil to slide this switch toward the middle of the sensor to switch it ON.
- 3 As soon as you activate the sensor via the switch (ON), seizure detection will begin. Two green PPG LEDs will light up brightly on the dark underside. One LED may turn off automatically when there is too much environmental light.
- 4 You can now connect the sensor to the charger and plug it into a power supply to stop seizure detection.

Disconnecting the sensor from the charger starts the seizure detection.

If you wish to turn the sensor OFF to prevent the battery from running out, use the ON/OFF switch to turn the sensor OFF. The sensor cannot charge when it is switched OFF.

Step 3: Using the alarm station


Placement

Position the alarm station where you want to receive alarms. It is important that the sound signals from the alarm station are clearly audible. The recommended location is in the bedroom of the operator or caregiver.

The distance between the alarm station and a worn sensor is limited and depends on the structure of the building where NightWatch+ is used. The typical range is approximately 15 meters. Be aware that wearing the sensor and especially blocking the signal with your body decreases the range. The alarm station will trigger a technical alarm if the sensor is out of range. If this happens, move the alarm station closer to the sensor, or ensure that the signal does not have to pass through too many walls.


Start-up

Ensure that the activated sensor is connected to the charger and that the alarm station is plugged into a power supply. Approximately 30 seconds after connecting the alarm station to a power supply, blue LED 3 (●) on the alarm station will begin blinking ('sensor charging'). Once the sensor is fully charged, blue LED 3 will remain steadily lit.

Adjusting the volume

When no alarm is active, the volume of the alarm station can be adjusted as needed using the volume controls (🔊) and (🔇) (see Section 5.3). Ensure that the sound signal from the alarm station is audible and loud enough to wake you while sleeping.

Your system is now ready to use.

Testing the alarm system

Remove the sensor from the charger and verify that within five seconds the alarm station begins blinking the technical alarm LED (orange) and emits an alarm sound indicating that no heart rate has been tracked yet. Notice that the alarm station responds when the sensor is removed from the charger, confirming that they are connected. When the sensor is reconnected to the charger, the alarm will stop within five seconds.

CAUTION: Only use the power adapters as supplied by LivAssured. Using a different charger or cable could damage the device and/or affect its performance.

PLEASE NOTE: The alarm station power supply plug and the sensor charging supply plug are the disconnecting devices from the electricity from the outlet. Make sure that the power supplies are always accessible.

5.2 Daily use

1

Disconnect the sensor from the charger. NightWatch+ is now in standard detection mode, where rapid shaking movements or a low heart rate are required to trigger an epilepsy alarm.

2


The alarm station will begin blinking the technical alarm LED (orange) and emit an alarm sound indicating that no heart rate has been tracked yet. This signal also confirms that the sensor and alarm station are communicating properly and that the alarms are audible.

WARNING: If the alarm station does not emit an alarm when disconnecting the sensor from the charger, it will not be able to notify you when a seizure event is detected. Please contact LivAssured in this case. See section 12 for contact details.


3

Place the sensor directly on the skin around the upper arm, on the thickest part of the biceps, with the sensor facing forward to prevent the wearer from lying on it when sleeping on their side. Lying on the sensor could disrupt tracking the heart rate. It does not matter if the indicator light is facing up or down.


4

The alarm station and sensor will blink green when the heart rate is tracked.

5

The wearer will now go to sleep. Once the wearer is lying down and the sensor has detected very little to no movement for two minutes, the more sensitive epilepsy detection algorithms become active. When this occurs, the blinking green LED changes to a steady green LED. NightWatch+ is now in rest mode and will also produce epilepsy alarms when vibrating movements or heart rate increases occur.


As soon as the wearer rises or gets out of bed, NightWatch+ automatically switches back to standard mode until the wearer lies down with no movement again for two minutes.

6

When the wearer wakes up, place the sensor back on the charger to stop detection. Once the sensor is fully charged, it will be ready for a new daily use cycle.

See chapter 6 for a list of all NightWatch+ alarms

PLEASE NOTE: When the sensor is switched ON and removed from the charger, both green PPG LEDs will light up. Both green LEDs should have the same light intensity. After a while, one of the green LEDs may switch itself off due to environmental light to optimize the heart rate tracking, this is normal.

PLEASE NOTE: As soon as the sensor is disconnected from the charger, the alarm station will emit an alarm sound until the heart rate has been tracked. Once the heart rate has been tracked, LED 4 (green) will start to blink or light up continuously. Check this regularly.

WARNING: Please ensure the alarm station is not covered during use as this can result in inaudible alarms.

WARNING: Damage to or degradation of the PPG sensor can result in the sensor not tracking heart rate. If the sensor cannot track the heart rate while wearing the sensor, please contact LivAssured. See section 12 for contact details.


WARNING: Please ensure that the alarm station is powered before the sensor is removed from the charger or you won't be able to hear if the alarm station is properly working.

5.3 Managing alarm sounds

The volume buttons and allow you to adjust the volume of the alarms emitted by the alarm station as well as to pause the audio of an alarm. You can change the volume when no alarm is currently active.

Adjusting the volume

- The buttons and can be used to adjust the alarm volume. You will hear the epilepsy alarm volume increase or decrease. Release the button at the desired volume. The selected volume setting is now saved and will remain unchanged even if the alarm station is disconnected from the power supply.

Pausing the audio of an alarm

- Either button can be pressed to pause the audio during an alarm.
- When the alarm audio is paused, a flashing green LED will appear on the alarm station.
- For epilepsy alarms, the flashing green light indicates that you have acknowledged the epilepsy alarm. Seizure detection will automatically resume three minutes after the onset of the alarm.

Resuming the audio of an alarm

- You can manually resume a paused alarm by pressing either button again.
- If a technical alarm is paused but remains active after 10 minutes, the alarm audio will automatically resume. For example, if you pause a technical alarm due to 'no heart rate found', the alarm will sound again if no heart rate is detected after 10 minutes.

WARNING: Ensure that the alarm sound of the alarm station is set to an audible sound level during use, or you will not be notified of an alarm sound.

CAUTION: Very loud sounds can trigger epileptic seizures in some people. Make sure that the signals can be heard clearly by the person who needs to respond to them but are not too loud for the person wearing the sensor.

5.4 Charging the sensor

The sensor is powered by a rechargeable battery and can be charged by connecting the supplied USB-C power adapter to the charging port. When the sensor is connected to the charger, the green PPG LEDs will turn off. The sensor is designed to remain on the charger for extended periods. When NightWatch+ is not in use, leave the sensor switched ON and connected to the charger. It will automatically stop charging when the battery is full and cannot overcharge. Seizure detection is disabled during charging.

The sensor will only charge when switched ON. The sensor will not charge if the ON/OFF switch is set to OFF.

When the sensor is not in use and cannot be connected to a charger, for example during transport, the ON/OFF switch should be set to OFF to prevent battery depletion.

Sensor battery life:

- ▼ Battery charging time: Approximately 2 hours.
- ▼ Fully charged battery life: At least 12 hours.

WARNING: Do not wear the sensor while it is charging because seizure detection is disabled while charging.

CAUTION: Contact LivAssured to replace the sensor if a fully charged battery of the sensor is empty before the end of a single cycle of use.

CAUTION: Always keep the sensor connected to the charger until use or the device may not be able to detect seizures during the full use cycle.

CAUTION: The sensor cannot charge when it is switched OFF. Always leave the sensor switched ON when it is connected to the charger.

PLEASE NOTE: The system starts automatically as soon as the sensor is removed from the charger. It is not possible to 'overcharge' the sensor.

PLEASE NOTE: Charging the sensor in the maximum operating temperature of 35°C may cause the contact surface to become 45°C, do not hold the sensor for longer than 10 minutes if it feels hot. Keep the sensor away from radiators and other sources of heat.

5.5 Maintenance & cleaning

NightWatch+ does not require periodic maintenance. However, for hygiene reasons, the sensor should be cleaned regularly with a damp cloth and disinfectant. Replace the comfort patch when it becomes dirty, damaged, or starts to come loose.

CAUTION: Clean the sensor regularly to reduce the risk of the wearer developing skin irritation.

CAUTION: Do not use excessive water to clean the device. Do not submerge the device or parts of the device in water. Do not wear the device in the rain, bath or in shower. This can damage the device.

5.6 Reuse

If NightWatch+ is to be used by a different person, the sensor should be cleaned with a damp cloth and disinfectant. A new piece of elastic strap should be used to fit the sensor to the new user, following the instructions in Section 0.

Data from the device is stored only in the NightWatch Portal if the device was connected to the internet. The data stored in the Portal consists of recorded device data and is not linked to any individual. If the Portal was used, you can delete the data before the device will be used by another person. See Chapter 7.3 for instructions on how to delete data from the Portal.

CAUTION: If the sensor is to be worn by a different user, it is recommended to be cleaned with a damp cloth and a disinfectant or else skin irritation can occur.

5.7 Transport or storage

When storing NightWatch+ or taking it along while traveling, both the alarm station and sensor should be switched OFF to prevent battery depletion. Disconnect the power supply from the alarm station. The "power lost" alarm will sound. Press one of the volume buttons to stop the 'power lost' alarm, and the alarm station will turn off. On the sensor, use a pointed object or a pen/pencil to slide the ON/OFF switch to OFF to prevent battery depletion.

CAUTION: The sensor and alarm station should be powered OFF during transport or storage or else the battery may be damaged. If you wish to use the system again, you can use the ON/OFF switch to switch the sensor back ON and charge it.

6 Signals from the sensor and the alarm station

Both the sensor and alarm station display signals indicating the system's status. The sensor has a light in the notch on the (white) top. The brightness of this indicator is optimized for use in a darkened room and may be more difficult to see in daylight or a brightly lit area. The alarm station uses both light and sound signals. These signals are explained in the following section.

-  LED 1 (blue) – Power
-  LED 2 (blue) – Network connection
-  LED 3 (blue) – Sensor charging status
-  LED 4 (green) – Audio paused state / Operation mode
-  LED 5 (orange) – Technical alarm
-  LED 6 (red) – Epilepsy alarm

6.1 Alarm signals

Light Indicator	Pausing alarm audio signals
	Flashing Except for the “no power” alarm, all alarm audio signals can be paused by pressing one of the volume buttons on the alarm station. A flashing green light, along with the existing alarm lights, will indicate that the alarm is paused. New alarms will be audible. <i>Example: Pressing one of the volume buttons while the “out of range” alarm is sounding will cause the alarm station lights to display as shown above. The flashing green audio-paused signal light will indicate that the alarm is paused, while the orange alarm light will continue to blink twice repeatedly until the issue is resolved.</i>

Light indicator	Sound	Meaning
	Red Flashing Fast beeping melody 	Epileptic seizure detected Check on the wearer and provide assistance as needed, following your doctor's instructions. Press either button to acknowledge the alarm and pause the audio. Seizure detection will automatically resume after three minutes.
	Orange Blinking 1x Monotone melody of three beeps 	Unable to track heart rate Check if the sensor is worn correctly. Press either button to pause the audio.
	Orange Blinking 2x Falling melody of three beeps 	Out of range The distance between the alarm station and the sensor is too great, or the sensor is switched OFF. Press either button to pause the audio.
	Blue Blinking 4x Orange Blinking 4x Rising melody of three beeps 	Sensor battery low Charge the sensor. Press either button to pause the audio.
	All off Continuous tone 	The alarm station has no power Check the power connection. Note: Pressing either button will turn the alarm station OFF.

6.2 Information signals

Light indicator	Sound	Meaning
 Green Blinking 	None	The heart rate is tracked
 Green Continuously on 	None	The heart rate is tracked while the wearer has been at rest The wearer remained at rest (no movement) in a horizontal position for two minutes. NightWatch+ is now operating in rest mode for seizure detection.
 Blue/Green Blinking 	None	The sensor is connected to the charger and is charging
 Blue/Green Continuously on 	None	The sensor is connected to the charger and is fully charged. Disconnect the sensor from the charger to begin seizure detection.

Operating modes

NightWatch+ has two modes for seizure detection: standard mode and rest mode.

Standard mode

When the sensor is removed from the charger, NightWatch+ operates in standard mode for seizure detection. Seizure alarms will be triggered by shaking movements or a detected low heart rate.

Rest mode

When the wearer's heart rate is tracked in a horizontal position (angle lower than 45 degrees) and minimal or no movement is detected for two minutes, the green indicator light will remain continuously on. This indicates that NightWatch+ is operating in rest mode.

In rest mode, seizure detection becomes more sensitive. NightWatch+ will also respond to smaller movements and increases in heart rate. If the sensor detects a vertical position, for example when the wearer gets out of bed, or when an alarm sounds, NightWatch+ will switch back to standard mode.

CAUTION: Be aware that if the sensor is not horizontal during sleep, NightWatch+ will stay in standard mode increasing the possibility that the caregiver is not warned of an epileptic seizure while using NightWatch+.

6.3 Responding to alarms

Simultaneous alarm conditions

NightWatch+ epilepsy and technical alarms are assigned a priority. If multiple alarm conditions occur at the same time, NightWatch+ announces the highest priority alarm.

Epilepsy alarm


● **Blinking red: Epilepsy alarm**

♪ **Sound: Fast beeping melody**

Priority: high

Epilepsy alarm!

A possible epileptic seizure has been detected. Check on the wearer and, if necessary, provide aid as instructed by your physician. This alarm will remain active until someone presses one of the volume buttons. When pressed, NightWatch+ will automatically resume seizure detection within a maximum of three minutes.


NightWatch+ may have triggered an epilepsy alarm for one of the following reasons:

Type	Description
Rhythmic movements	Rhythmic movements have been detected, which could indicate epileptic activity.
Vibrating movements	Vibrating or muscle tension movements have been detected, which could indicate epileptic activity.
Shaking movements	Intense shaking movements have been detected, which could indicate epileptic activity.
Heart rate increase	The average tracked heart rate rose suddenly, which could indicate epileptic activity.
High heart rate	The tracked heart rate is significantly higher than before, which could indicate epileptic activity.
Low heart rate	The tracked heart rate is low, which could indicate epileptic activity.

Technical alarms

● **Blinking orange: Technical alarm**

There are 3 types of technical alarms:

Sound	Alarm LED 5 (orange)	Priority	Technical alarm
1 Falling melody 	Blinking 2x repeatedly	High	Sensor out of range or sensor switched off
2 Monotone melody 	Blinking 1x repeatedly	Medium	Sensor unable to track heart rate
3 Rising melody 	Blinking 4x repeatedly	Low	Sensor battery low

PLEASE NOTE: If the alarm station makes a sound not described in the manual, please contact LivAssured. See section 12 for contact details.

1. Out of range

There is no connection between the sensor and the alarm station. Several causes are possible:

- ▼ The sensor and alarm station are too far apart. Move the alarm station and sensor closer together.
- ▼ The sensor is switched OFF. Turn the sensor ON by sliding the ON/OFF switch. At least one green LED at the bottom of the sensor should illuminate when the sensor is ON and disconnected from the charger.
- ▼ The alarm station is not connected to a power supply.
- ▼ The sensor displays a continuous (non-blinking) orange LED in the notch. If this occurs, press and hold both volume buttons on the alarm station for 20 seconds. The sensor should automatically reconnect to the alarm station.

2. Unable to track heart rate

You will hear this alarm as soon as you remove the sensor from the charger. The alarm will stop as soon as NightWatch+ tracks a heart rate. Hold still to allow the heart rate to be tracked more quickly. This alarm will sound again if the sensor fails to track both heart rate and movement for at least two minutes. Please note that NightWatch+ is not a heart rate monitor but uses heart rate tracking alongside movement detection to notify for seizures.

There are several possible causes of a temporary or permanent failure to track a reliable heart rate:

- ▼ The PPG sensor is not correctly positioned on the skin. Ensure that the green LEDs on the underside of the sensor make direct contact with the skin.
- ▼ The sensor may be too loose, have slid off, or has been removed. If so, tighten the elastic strap slightly.
- ▼ The wearer may be lying on the sensor, making heart rate detection difficult. Adjust the sensor position to prevent the wearer from lying on top of it. Try placing it on top of the bicep muscle, closer to the inside of the arm (near the armpit) rather than the outer side.

3. Sensor battery low

This alarm indicates that the sensor battery is nearly empty. NightWatch+ will no longer function until the sensor is charged. Connect the sensor to the charger.


If this alarm persists despite correct and sufficient charging, please contact 'Service and support'. See Section 12 for contact details.

Power loss alarm

- No LED on
- ♪ Sound: constant beep: Power loss

Priority: high

When the alarm station loses power, a 'power loss' alarm will sound. Restore power to the alarm station to stop the alarm, or press one of the volume buttons to turn off the alarm station completely.

7 Seizure monitoring

7.1 The NightWatch Portal

NightWatch+ does not need to be connected to the internet to function correctly. However, if it is connected to the internet during use, the system will transmit captured data and alarms to the NightWatch Portal. This data can be accessed via portal.nightwatchepilepsy.com.

The Portal provides insight into events that occurred during sleep, helping users better understand why NightWatch+ did or did not trigger alarms. Additionally, it can be used to view alarm frequency and share the device data, along with personal notes, with a neurologist. It is also useful for troubleshooting.

The alarm station does not store any data internally. It will only send data to the Portal when it is actively recording and connected to the internet via a wired connection.


7.2 How to connect to the NightWatch Portal

To connect to the NightWatch Portal, the NightWatch+ alarm station must be permanently wired to an internet access point, such as a router.

- 1 Locate the nearest internet access point (router) to the NightWatch+ alarm station and connect them using an Ethernet cable.


- 2 If the nearest access point is too far for a direct cable connection, we recommend using a powerline adapter, or mobile router to bridge the connection.

- 3 Once NightWatch+ is connected to the internet, the second blue indicator LED on the alarm station will illuminate, confirming the internet connection.

- 4 After a recording has been made while NightWatch+ was connected to the internet, your device can be added to your Portal account. Create your NightWatch Portal account at portal.nightwatchepilepsy.com and follow the steps to add your device.

- 5 Additional information about the Portal and its features can be found by scanning the QR code or visiting nightwatchepilepsy.com/helpdesk

7.3 Does the NightWatch Portal respect my privacy?

If your NightWatch+ remains connected to the internet during use, it will automatically send recordings to the NightWatch Portal. This data includes information about heart rate, movements, and alarms. However, it does not contain any personally identifiable information about the wearer and is stored encrypted and anonymously. NightWatch Portal does not require users to enter personal data such as names or addresses. You may request our customer support team to delete your device data at any time, or you can remove your own historic recordings from the settings page of the Portal.

7.4 Connecting NightWatch+ to alarm systems

NightWatch+ offers two different methods for connecting to information or alarm systems. These options are available only for professional use upon request.

- The RJ-11 port allows connection to a Distributed Information System
- The RJ-45 port allows connection to a Distributed Alarming System.

No other interconnections are possible. For more information, please contact LivAssured. See Section 12 for contact details.

8 Specifications

8.1 Technical Specifications

Operating mode	[Sensor]: Body worn
Weight	[Sensor]: 35g [Alarm station]: 90g
Dimensions (L x W x H)	[Sensor]: 72 mm x 52 mm x 14 mm [Alarm station]: 100 mm x 100 mm x 28 mm
Supply voltage	[Sensor]: FRIWO NEO006.01-X-05, 100V–240V AC, 50Hz–60Hz, IEC60601-1 protection class II, 5VDC/1.4A. Lifetime: 20 years at continuous use. [Alarm station]: FRIWO NEO006.01-X-05, 100V–240V AC, 50Hz–60Hz, IEC60601-1 protection class II, 5VDC/1.4A. Lifetime: 20 years at continuous use.
Current consumption	[Sensor]: 0.1A (RMS) max. [Alarm station]: 0.2A (RMS) max.
Internal transmitters	DECT, operating frequency (send/receive) 1880–1900MHz, 23 dBm
Casing protection	[Sensor]: IP21 – This means that the device is protected against solid foreign objects of 12,5 mm diameter and greater, and against vertically falling water drops. [Alarm station]: IP22 – This means that the device is protected against solid foreign objects of 12,5 mm diameter and greater, and against vertically falling water drops when enclosure tilted up to 15°
Applied part	Sensor, type BF
Battery	[Sensor] Built-in Lithium-ion battery 3.7V, 450mAh, not replaceable [Alarm station] Built-in Lithium-ion battery 3.7V, 450mAh, not replaceable Both batteries comply to IEC62133-2:2017 and UN38.3
Sound pressure range	40 to 80 dBA for all alarms (1m radius)
Alarm condition delay	2 seconds max.
Pulse rate accuracy	(30 – 210 bpm) ± 1.31 bpm (RMS)
Pulse rate acc. method	Electronic pulse simulator

CAUTION: Batteries inside NightWatch+ cannot be replaced. Trying to replace the batteries can damage NightWatch+ which could result in its incorrect operation. Under normal use conditions the batteries service life is 5 years. Please contact LivAssured if you think the battery is not working properly. See section 12 for contact details.

8.2 Environmental conditions

Operating conditions

- ▼ Temperature range of +5°C to +35°C
- ▼ Relative humidity range of 15% to 90%, non-condensing, but not requiring a water vapor partial pressure greater than 50hPa
- ▼ Atmospheric pressure range of 700hPa to 1060hPa

PLEASE NOTE: Charging the sensor in the maximum operating temperature of 35°C may cause the contact surface to become 45°C, do not hold the sensor for longer than 10 minutes if it feels hot. Keep the sensor away from radiators and other sources of heat.

CAUTION: Using NightWatch+ in an environment above 35°C may cause the contact surface of the sensor to become hot and unsuitable to wear.

Transport and storage condition limits

- ▼ Temperature range of -25°C to +70°C
- ▼ Relative humidity range of 15% to 90%, non-condensing
- ▼ Atmospheric pressure range of 700hPa to 1060hPa

Recommended storage conditions

- ▼ Temperature range -20°C to +25°C for a maximum of 3 months with the sensor turned off

CAUTION: Storage temperatures above 25°C will increase the rate of self-discharge, reducing the available capacity of the battery. A reduced capacity of the battery may cause the sensor to be empty before the end of a single cycle of use and will not be able to detect seizures during the full use cycle.

8.3 Electromagnetic environment conditions

NightWatch+ is intended for use in the electromagnetic environment specified below. The operator should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The above listed models use RF energy only for its internal function. Therefore, its RF emissions are ultra low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The above listed models are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Enclosure Port			
Immunity test	Test Condition	IEC 60601 Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV Contact ± 2,4,8,15 kV Air	±8 kV Contact ± 15 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF EM fields and Proximity fields from RF wireless communications equipment IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM 1kHz	10 V/m 80 MHz – 2,7 GHz	Mains power quality should be that of a professional healthcare facility environment and home healthcare environment.
	385MHz (18Hz Pulse Modulation)	27 V/m	
	450MHz (FM+/-5KHz deviation 1kHz sine or 18Hz Pulse Modulation)	28 V/m	
	710MHz (217Hz PM)	9 V/m	
	745MHz (217Hz PM)	9 V/m	
	780MHz (217Hz PM)	9 V/m	
	810MHz (18Hz PM)	28 V/m	
	870MHz (18Hz PM)	28 V/m	

	930MHz (18Hz PM)	28 V/m	
	1720MHz (217Hz PM)	28 V/m	
	1845MHz (217Hz PM)	28 V/m	
	1970MHz (217Hz PM)	28 V/m	
	2450MHz (217Hz PM)	28 V/m	
	5240MHz (217Hz PM)	9 V/m	
	5500MHz (217Hz PM)	9 V/m	
	5785MHz (217Hz PM)	9 V/m	
RATED power frequency magnetic fields IEC 61000-4-8	50Hz or 60Hz	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Proximity magnetic fields IEC 61000-4-39	30 kHz (CW)	8 A/m	Proximity magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
	134,2 kHz (2,1 kHz PM)	65 A/m	
	13,56 MHz (50 kHz PM)	7,5 A/m	

Input a.c. power PORT			
Immunity test	Test Condition	IEC 60601 Compliance level	Electromagnetic environment - guidance
Electrical fast transient/bursts IEC 610004-4	± 2 kV 100kHz Repetition frequency	± 2 kV	Mains power quality should be that of a professional healthcare facility environment and home healthcare environment.
Surges IEC 61000-4-5	± 0,5kV, ±1 kV line(s) to line(s)	± 1 kV, Differential mode	Mains power quality should be that of a professional healthcare facility environment and home healthcare environment.
Conducted RF induced by RF fields IEC 61000-4-6	3 Vrms 150 kHz - 80 MHz also 6 Vrms ISM and Amateur Radio Bands a) 80% AM 1kHz	3 Vrms 150 kHz - 80 MHz also 6 Vrms ISM Radio Bands a) 80% AM 1kHz	Mains power quality should be that of a professional healthcare facility environment and home healthcare environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U _r ; 0°,45°,90°,135°,180°,225°,270°,315° 0% U _r ; 0° 0% U _r ; 70% 0% U _r ; 0%	0,5 Cycles 1Cycle 25/30 Cycles (50/60Hz) 250/300 Cycles (50/60Hz) (5s)	Mains power quality should be that of a professional healthcare facility environment and home healthcare environment. If the user of the above listed models requires continued operation during power mains interruptions, it is recommended that the

			above listed models are powered from an uninterruptible power supply or battery.
Comment: a) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.			

Signal input/output parts PORT			
Immunity test	Test Condition	IEC 60601 Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV Contact ± 2,4,8,15 kV Air	±8 kV Contact ± 15 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/bursts IEC 610004-4	± 1 kV 100kHz Repetition frequency	± 1 kV	Mains power quality should be that of a professional healthcare facility environment and Home healthcare environment.

Conducted RF induced by RF fields IEC 61000-4-6	3 Vrms 150 kHz - 80 MHz also 6 Vrms ISM and Amateur Radio Bands a) 80% AM 1kHz	3 Vrms 150 kHz - 80 MHz also 6 Vrms ISM and Amateur Radio Bands a) 80% AM 1kHz	Mains power quality should be that of a professional healthcare facility environment and Home healthcare environment.
Comment: a) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.			

CAUTION: Using NightWatch+ in electromagnetic environments outside these descriptions may cause disturbances leading to a loss of connection between the sensor and the alarm station.

8.4 Regulatory and Compliance

NightWatch+ is in compliance with the following standards for medical devices and radio equipment and has been subject to specific laboratory testing to assess its safety, electromagnetic compatibility, usability, and biocompatibility.

Basic safety and essential performance	EN IEC 60601-1:2005+A1:2012+A2:2020 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
EMC requirements	EN IEC 60601-1-2:2014+A1:2020 – Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests EN 301 489-1 V2.2.3 ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements; Harmonised Standard for ElectroMagnetic Compatibility EN 301 489-6 V2.2.1 ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 6: Specific conditions for Digital Enhanced Cordless Telecommunications (DECT) equipment;
DECT	ETSI EN 301 406 V2.2.2 (2016-09) Digital Enhanced Cordless Telecommunications (DECT); Harmonised Standard for access to radio spectrum; Part 1: DECT, DECT Evolution and DECT ULE
Home care use	EN IEC 60601-1-11:2015+A1:2020 Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance – Collateral standard: Requirements for medical electrical equipment and medical electrical Systems used in the home healthcare environment
Medical Alarm systems	EN IEC 60601-1-8:2006+A1:2012+A2:2020 Medical Electrical equipment Part 1-8: General requirements for basic safety and essential performance – Collateral standard: General Requirements, tests and guidance for alarm systems in medical electrical equipment.
Biocompatibility	EN ISO 10993-1:2020 Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process

EN ISO 10993-5:2009 Biological evaluation of medical devices, Part 5: EN Tests for in vitro cytotoxicity
EN ISO 10993-10:2023 Biological evaluation of medical devices, Part 10: Tests for skin sensitization
EN ISO 10993-23:2021 Biological evaluation of medical devices, Part 23: Tests for irritation

Risk Management	EN ISO 14971:2019+A11:2021 Application of risk management to medical devices
Medical Device Software	EN IEC 62304:2006+A1:2015 Software life-cycle processes
Cybersecurity	EN IEC 81001-5-1:2021 Health software and health IT systems safety effectiveness and security Part 5-1: Security Activities in the product lifecycle AAMI TIR57:2016/(R)2023 Principles for medical device security – Risk management
Labelling and symbols	EN ISO 15223-1:2021 Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirement
Usability	EN IEC 60601-1-6:2010+A1:2013+A2:2020 Medical electrical equipment Part 1-6 General requirements for basic safety and essential performance– Collateral standard: Usability EN 62366-1:2015+A1:2020, Medical devices -- Part 1: Application of usability engineering to medical devices

8.5 Cybersecurity

LivAssured works according to the latest standards for security risk management to keep your data safe. The alarm station and sensor do not store any health data internally. Data recorded by the sensor is recorded for no more than 5 minutes.

The alarm station and the sensor are connected via DECT Ultra Low Energy (ULE). DECT has been standardized for the purpose of cordless telephony and ULE provides enhanced encryption. A connection according to the DECT standard uses subscription, authentication and encryption techniques to secure the data stream

and uses a dedicated radio frequency for high stability. The DECT connection used by the NightWatch+ has been tested according to the latest standards for DECT.

Connecting to the NightWatch Portal is entirely optional, the system is fully functional without this connection. When connected, communication between the alarm station and the Portal is encrypted. The data is stored at a partner with ISO 27001:2022 certification, an international security standard that ensures having and maintaining a high security level of the data. Moreover, the data that the Portal displays is anonymous device recorded data and cannot be identified to any person.

9 Incident reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the country in which the user is established.

10 Service life and guarantee

NightWatch+ has a 2-year guarantee. In the event NightWatch+ is not working or seems to be working incorrectly, please contact LivAssured. See section 12 for contact details. The expected service life of NightWatch+ including the parts shipped with the device in case of daily use is 5 years.

11 Disposal

At the end of its useful life, NightWatch+ (with its battery) must be disposed of in accordance with local law and the local code concerning the disposal of electrical and electronic equipment including lithium-ion batteries. Do not discard NightWatch+ in a standard trash bin.

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Support: +3185 0601252

Email: info@nightwatchepilepsy.com

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